A SWISS BRAND

**BLOOD
GLUCOSE
MONITOR
SYSTEM**

**MODEL MINI
USER'S MANUAL**

www.dr-frei.com

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DEAR DR. FREI MINI SYSTEM OWNER:

Thank you for your purchase of blood glucose monitoring system **Dr. Frei MINI**. We are sure that having appraised worthily the high quality and reliability of this device you will become a regular user of the products of Swiss Trademark **Dr. Frei®**.

This manual provides important information to help you to use the system properly. Before using this product, please read the manual through carefully.

Regular monitoring of your blood glucose levels can help you and your doctor to gain better control of your diabetes. Due to its compact size and easy operation, you can use the **Dr. Frei MINI** Blood Glucose Monitoring System to easily monitor your blood glucose levels by yourself anywhere and anytime.

For all questions concerning the device please contact your local distributor or **Dr. Frei®** service center in your country.

TABLE OF CONTENTS

| | |
|---|-----------|
| SAFETY REGULATIONS. READ BEFORE TESTING..... | 6 |
| Important Information | 7 |
| Intended Use | 7 |
| Test Principle | 7 |
| Contents of System | 8 |
| Meter Overview | 9 |
| Display Screen | 10 |
| Test Strip | 11 |
| Setting the Meter | 12 |
| TESTING WITH BLOOD SAMPLE..... | 14 |
| Preparing the Lancing Device for Blood Testing | 14 |
| Preparing the Puncture Site | 14 |
| Alternative Site Testing | 16 |
| Performing a Blood Glucose Test..... | 17 |
| METER MEMORY | 19 |
| Reviewing Test Results | 19 |
| Result Readings | 20 |
| Comparing Meter and Laboratory Results..... | 21 |
| MAINTENANCE | 22 |
| Battery | 22 |
| Caring for Your Meter | 24 |
| Caring for Your Test Strips..... | 24 |
| Error Messages | 24 |
| INFORMATION FOR THE SERVICE CENTRE EXPERT..... | 26 |
| Important Control Solution Information..... | 27 |
| Performing a Control Solution Test..... | 28 |
| Troubleshooting..... | 30 |
| SYMBOL INFORMATION | 32 |
| SPECIFICATIONS..... | 33 |

IMPORTANT SAFETY INSTRUCTIONS.READ BEFORE USE!

1. Use this device **ONLY** for the intended use described in this manual.
2. Do **NOT** use accessories which are not specified by the manufacturer.
3. Do **NOT** use the device if it is not working properly or if it is damaged.
4. Do **NOT** use the equipment in places where aerosol sprays are being used or where oxygen is supplying.
5. Do **NOT** under any circumstances use the device on newborns or infants.
6. This device does **NOT** serve as a cure for any symptoms or diseases. The data measured is for reference only.
7. Before using this device to test blood glucose, read all the instructions thoroughly and practice the test. Carry out all the quality control checks as directed.
8. Keep the device and testing equipment away from young children. Small items such as the battery cover, batteries, test strips, lancets and vial caps are choking hazards.
9. The usage of this instrument in a dry environment, especially if synthetic materials are present (synthetic clothing, carpets etc.) may provoke damaging static discharges that cause erroneous results.
10. Do not use this instrument in close proximity to sources of strong electromagnetic radiation, as these may interfere with the accurate operation.

**KEEP THESE INSTRUCTIONS IN A SAFE PLACE
FOR FUTURE REFERENCE**

Important Information

- Severe dehydration and excessive water loss may cause readings which are lower than actual values. If you believe you are suffering from severe dehydration, consult a healthcare professional immediately.
- If your blood glucose results are lower or higher than usual, and you do not have any symptoms of illness, first repeat the test. If you have symptoms or continue to get results which are higher or lower than usual, follow the treatment advice of your healthcare professional.
- Use only fresh whole blood samples to test your blood glucose.
- Using other substances will lead to incorrect results.
- If you are experiencing symptoms that are inconsistent with your blood glucose test results and you have followed all the instructions given in this owner's manual, contact your healthcare professional.
- We do not recommend using this product on severely hypotensive individuals or patients in shock. Readings which are lower than actual values may occur for individuals experiencing a hyperglycaemic-hyperosmolar state, with or without ketosis.
- Please consult the healthcare professional before use.

Intended Use


This system is intended for use outside the body (in vitro diagnostic use) by people with diabetes at home and by health care professionals in clinical settings as an aid to monitoring the effectiveness of diabetes control. It is intended to be used for the quantitative measurement of glucose (sugar) in fresh whole blood samples (from the finger, palm, forearm, upper arm, calf and thigh). It should not be used for the diagnosis of diabetes, or testing on newborns.

Test Principle

Your system measures the amount of sugar (glucose) in whole blood. The glucose testing is based on the measurement of electrical current generated by the reaction of glucose with the reagent of the strip. The meter measures the current, calculates the blood glucose level, and displays the result. The strength of the current produced by the reaction depends on the amount of glucose in the blood sample.

CONTENTS OF SYSTEM

Your new **Dr. Frei MINI** system kit includes:


- ① Meter.
- ② Owner's Manual.
- ③ Storage Case.
- ④ Warranty Card.
- ⑤ Lancing Device with One Clear Cap.
- ⑥ Lithium battery CR 2032 3V.
- ⑦ Test Strips.
- ⑧ Sterile Lancets.

NOTE

All items can be purchased separately and some accessories may not be included in the kit. If you wish to purchase any accessories, please contact our local customer service.

METER OVERVIEW

① **Test Slot**

Insert test strip here to turn the meter on for testing.

② **Display Screen**

③ **M Button**


Enter the meter memory and silence a reminder alarm.

④ **Battery**

⑤ **SET Button**



Enter and confirm the meter settings.

DISPLAY SCREEN

| | |
|---|---|
|  | Battery indicator icon |
|  | Symbol that shows that you may apply a drop of blood on the test-strip |
|  | Test Result (is gotten after configuration procedure and the test itself) |
|  | Symbol of test-results which are keeping in the memory of device. |
|  | Time |
| mmol/L | Measurement Unit |
|  | Error Message |

TEST STRIP

ATTENTION:

The front side of test strip should face up when inserting test strip.
Test results might be wrong if the contact bar is not fully inserted into the test slot.

NOTE


The Dr. Frei MINI meter should only be used with Dr. Frei MS test strips.
Using other test strips with this meter can produce inaccurate results.

SETTING THE METER

Before using your meter for the first time or if you change the meter battery, you should check and update these settings. Make sure you complete the steps below and have your desired settings saved.

Entering the Setting Mode.

Start with the meter off (no test strip inserted). Press SET to turn on the meter.

Step 1 Setting the date.

With the year flashing, press and hold the button **M** until the correct year appears. Press SET.

With the month flashing, press and hold the button **M** until the correct month appears. Press SET.

With the day flashing, press and hold the button **M** until the correct day appears. Press SET.

Step 2 Setting the time format.

Press and release **M** to select the desired time format - 12h or 24h. Press SET.

Step 3 Setting the time.

With the hour flashing, press and hold the button **M** until the correct hour appears. Press SET.


With the minute flashing, press and hold the button **M** until the correct minute appears. Press SET.

Step 4 Deleting the memory.

You may set up any or all of the reminder alarms (1-4). The meter displays "On" or "OFF" and "1:00", press **M** to turn on or turn off to set the first reminder alarm.

Press **M** to select "On", then press SET to set the hour. When the hour is flashing, press **M** to add an hour. Press SET to confirm and go to minutes, press **M** to add one minute. Hold **M** longer to add faster. Press SET to confirm and go to the next alarm setting.

Congratulations! You have completed all settings!

NOTE

- These parameters can ONLY be changed in the setting mode.
- If the meter is idle for 3 minutes during the setting mode, it will switch off automatically

TESTING WITH BLOOD SAMPLE

NOTE

To reduce the chance of infection:

- Never share a lancet or the lancing device.
- Always use a new, sterile lancet. Lancets are for single use only.
- Avoid getting hand lotion, oils, dirt, or debris in or on the lancets and the lancing device.

Preparing the Lancing Device for Blood Testing

Please follow the instructions in the lancing device insert for collecting a blood sample.


Preparing the Puncture Site

Please follow the suggestions below before obtaining a drop of blood:

- Wash and dry your hands before starting.
- Select the puncture site either at fingertips or another body parts (please see section “Alternative Site Testing” (AST) on how to select the appropriate sites).
- Clean the puncture site using cotton moistened with 70% alcohol and let it air dry.
- Rub the puncture site for about 20 seconds before penetration.
- Use a clear cap (included in the kit) while setting up the lancing device.


Fingertip Testing

Press the lancing device's tip firmly against the lower side of your fingertip. Press the release button to prick your finger, then a click indicates that the puncture is complete.

Blood from Sites Other Than the Fingertip

Replace the lancing device cap with the clear cap for AST. Pull the cocking control back until it clicks. When lancing the forearm, upper arm, hand, thigh, or calf, avoid lancing the areas with obvious veins because of excessive bleeding.

NOTE

- Choose a different spot each time you test. Repeated punctures at the same spot may cause soreness and calluses.
- Please consult your health care professional before you begin AST.
- It is recommended to discard the first drop of blood as it might contain tissue fluid, which may affect the test result.

Alternative Site Testing (AST)


Important: There are limitations with AST (Alternative Site Testing). Please consult your health care professional before you perform AST.

What is AST?

Alternative site testing (AST) means that people may use other than fingertips parts of the body to check their blood glucose levels. This system allows you to test on the palm, the forearm, the upper arm, the calf or the thigh with results equivalent to fingertip testing.

What is the advantage?

Fingertips feel pain more readily because they are full of nerve endings (receptors). At other body sites, since nerve endings are not so condensed, you will not feel as much pain as on the fingertips.

When to use AST?

Food, medication, illness, stress and exercise can affect blood glucose levels. Capillary blood at the fingertip reflects these changes faster than capillary blood at other sites. Thus, when testing blood glucose during or immediately after a meal, physical exercise, or any other event, **take a blood sample from your finger only.**

We strongly recommend that you perform AST ONLY at the following times:

- In a pre-meal or fasting state (more than 2 hours since the last meal).
- Two hours or more after taking insulin.
- Two hours or more after exercise.

Do NOT use AST if:

- You think your blood glucose is low.
- You are unaware of hypoglycemia
- You are testing for hyperglycemia
- Your AST results do not match the way you feel.
- Your routine glucose results are often fluctuate.

Performing a Blood Glucose Test

- Step 1** Wait for the meter to display the test strip "E" and blood drop "d" symbol.

- Step 2** Obtaining a blood sample.

Use the pre-set lancing device to puncture the desired site. After penetration, discard the first drop of blood with a clean cotton swab. Gently squeeze the punctured area to obtain another drop of blood.

Be careful NOT to smear the blood sample.


The volume of blood sample must be at least 0.5 microliter (μL) of volume.

- Step 3** Apply the sample.

Hold the blood drop to touch the absorbent hole of the test strip. Blood will be drawn in and after the confirmation window is completely filled, the meter begins counting down.


NOTE

- Do not press the punctured site against the test strip or try to smear the blood.
- If you do not apply a blood sample to the test strip within 3 minutes, the meter will automatically turn off. You must remove and reinsert the test strip to start a new test.
- The confirmation window should be filled with blood before the meter begins to count down. NEVER try to add more blood to the test strip after the drop of blood has moved away. Discard the used test strip and retest with a new one.
- If you have trouble filling the confirmation window, please contact your health care professional or the local customer service for assistance.

Step 4 Read your result.

The result of your blood glucose test will appear after the meter counts down to 0. This blood glucose result will automatically be stored in the memory.

Step 5 Eject the used test strip and remove the lancet.

To eject the test strip, point the strip at a disposal container for sharp items. The meter will switch itself off automatically after the test strip is ejected.

Step 6

Always follow the instructions in the lancing device insert when removing the lancet.

WARNING!

The used lancet and test strip may be biohazards.
Please discard them carefully according to your local regulations.

METER MEMORY

The meter stores the 20 most recent blood glucose test results along with respective dates and times in its memory. To enter the meter memory, start with the meter switched off.

Reviewing Test Results

Step 1 Press **M**.

"**M**" will appear on the display. Press **M** again, and the first reading you see is the last blood glucose result along with date, time and the measuring mode.


Step 2

Press **M** to recall the test results stored in the meter each time you press.

Step 3 Exit the meter memory.

After the last test results, press **M** again and the meter will be turned off.

NOTE

- Any time you wish to exit the memory, keep pressing **M** for 3 seconds or leave it without any action for 3 minutes. The meter will switch off automatically.
- Control solution results are NOT included in the day average.
- If using the meter for the first time, "----" displays when you recall the test results or review the average result. It indicates that there is no test result in the memory.

Result Readings

| APPEARS | WHEN GLUCOSE |
|---------|-----------------------------|
| Lo | < 20 mg/ dL (1.1 mmol/ L) |
| Hi | > 600 mg/ dL (33.3 mmol/ L) |

Comparing Meter and Laboratory Results

The meter provides you with plasma equivalent results. The result you obtain from your meter may differ somewhat from your laboratory result due to normal variation. Meter results may be affected by factors and conditions that do not affect laboratory results in the same way. To make an accurate comparison between meter and laboratory results, follow the guidelines below

- Take your meter with you to the lab.
- Wash your hands before obtaining a blood sample
- Use fresh capillary blood only.
- To compare the results use one and the same blood drop for testing.

You may still have a variation from the result because blood glucose levels can change significantly over short periods of time, especially if you have recently eaten, exercised, taken medication or experienced stress*. In addition, if you have eaten recently, the blood glucose level from a finger prick can be up to 70 mg/dL (3.9 mmol/L) higher than blood drawn from a vein (venous sample) used for a lab test**.

Factors such as the amount of red blood cells in the blood (a high or low hematocrit) or the loss of body fluid (dehydration) may also cause a meter result to be different from a laboratory result.

References:

*: Surwit, R.S., and Feinglos, M.N.: Diabetes Forecast (1988), April, 49-51.

** : Sacks, D.B.: "Carbohydrates." Burtis, C.A., and Ashwood, E.R.(ed.), Tietz Textbook of Clinical Chemistry. Philadelphia: W.B. Saunders Company (1994), 959.

MAINTENANCE

Battery

Your meter comes with lithium battery CR 2032 3V.

Low Battery Signal

If the meter is functional and the result remains accurate, but it is time to change the battery the “” symbol appears along with display messages.

Replacing the Battery

To replace the battery, make sure that the meter is turned off.

Step 1


Press the edge of the battery cover and lift it up to remove.

Step 2

Remove the old battery and replace with CR 2032 3V.

Step 3

Close the battery cover. If the battery are inserted correctly, you will hear a "beep" afterwards.

NOTE

- Replacing the battery does not affect the test results stored in the memory.
- As with all small batteries, these battery should be kept away from small children. If swallowed, promptly seek medical assistance.
- Battery might leak chemicals if unused for a long time. Remove the battery if you are not going to use the device for an extended period (i.e., 3 months or more).
- Properly dispose of the battery according to your local environmental regulations.

Caring for Your Meter

To avoid the meter and test strips attracting dirt, dust or other contaminants, please wash and dry your hands thoroughly before use.

Cleaning

- To clean the meter exterior, wipe it with a cloth moistened with tap water or a mild cleaning agent, then dry the device with a soft and dry cloth. Do **NOT** rinse with water.
- Do **NOT** use organic solvents to clean the meter.

Meter Storage

- Storage conditions: 4°C to 40°C (39.2°F to 104°F), below 95% relative humidity.
- Always store or transport the meter in its original storage case.
- Avoid dropping and heavy impact.
- Avoid direct sunlight and high humidity.

Caring for Your Test Strips

- Storage conditions: 4°C to 40°C (39.2°F to 104°F), below 85% relative humidity. Do **NOT** freeze.
- Store your test strips in their original vial only. Do not transfer to another container.
- Store test strip packages in a cool and dry place. Keep away from direct sunlight and heat.
- After removing a test strip from the vial, immediately close the vial cap tightly.
- Touch the test strip with clean and dry hands.
- Use each test strip immediately after removing it from the vial.
- Write the opening date on the vial label when you first opened it. Discard remaining test strips after 3 months.
- Do not use test strips beyond the expiry date. This may cause inaccurate results.
- Do not bend, cut, or alter a test strip in any way.
- Keep the strip vial away from children since the cap and the test strip may be a choking hazard. If swallowed, promptly see a doctor for help.

For further information, please refer to the test strip package insert.

Error Messages

| MESSAGE | WHAT DOES METER SAY? | WHAT IT MEANS | WHAT TO DO |
|---|--|--|--|
|  | Battery is dead. Please replace. | Appears when the batteries cannot provide enough power for a test. | Replace the batteries immediately. |
|  | Test strip has been removed during testing. | Remove the strip after applying blood to the absorbent hole. | Re-test with a new test strip. |
|  | The test strip has been used. Please replace. | Appears when a used test strip is inserted. | Repeat with a new test strip. |
|  | Room temperature is out of range. Thus unable to measure. | Appears when ambient temperature is below system operation range. | System operation range is 10°C to 40°C (50°F to 104°F). Repeat the test after the meter and test strip are in the above temperature range. |
|  | | Appears when ambient temperature is above system operation range. | |
|   | System error. | Problem with the meter. | Repeat the test with a new test strip. If the meter still does not work, please contact the customer service for assistance. |
|  | Memory error. | Problem with the meter. | Please contact your local customer service agent for help. |

INFORMATION FOR THE EXPERT OF THE LOCAL CUSTOMER SERVICE

CHECKING THE SYSTEM WITH CONTROL SOLUTION TM Dr. Frei

Attention! Control solution must be used only for 90 days from the day of its opening.

Before test starts:

Our control solution TM Dr. Frei contain a known amount of glucose that reacts with test strips. By comparing the result of your control solution test with the expected range printed on the test strip vial label, you can check that the meter and the test strips are working together as a system and that you are performing the test correctly. It is very important that you perform this simple check routinely to make sure you get accurate results. Control solution is not included in the system kit. You can test your system in our service centres indicated in warranty card.


Important Control Solution Information

- Use only Dr. Frei MINI control solutions with your meter.
- Do not use the control solution beyond the expiry date or 3 months after first opening. Write the opening date on the control solution vial and discard the remaining solution after 3 months.
- It is recommended that the control solution test be done at room temperature (20°C to 25°C / 68°F to 77°F). Make sure your control solution, meter, and test strips are at this specified temperature range before testing.
- Shake the vial before use, discard the first drop of control solution, and wipe off the dispenser tip to ensure a pure sample and an accurate result.
- Store the control solution tightly closed at temperatures between 2°C and 30°C (35.6°F and 86°F). Do NOT freeze.

Performing a Control Solution Test


Step 1 Insert the test strip to turn on the meter.

Insert the test strip into the meter.
Wait for the meter to display the test strip and blood drop symbols.

Step 2 Press **M** to mark this test as a control solution test.

If you press **M** again, the "C CL" will disappear and this test is no longer a control solution test.


Step 3 Apply control solution.

Shake the control solution vial thoroughly before use. Squeeze out a drop and wipe it off, then squeeze out another drop and place it on the tip of the vial cap.

Hold the meter to move the absorbent hole of the test strip to touch the drop. Once the confirmation window fills completely, the meter will begin counting down.

To avoid contaminating the control solution do not directly apply control solution onto a strip.

Step 4 Read and compare the result.

After counting down to 0, the control solution test result will appear on the display. Compare this result with the range printed on the test strip vial and it should fall within this range. If not, please read the instructions again and repeat the control solution test.

Results outside the range

If you continue to have test results fall outside the range printed on the test strip vial, the meter and strips may not be working properly.

WARNING!

The used lancet and test strip may be biohazards.
Please discard them carefully according to your local regulations

Troubleshooting

1. If the meter does not display a message after inserting a test strip:

| POSSIBLE CAUSE | WHAT TO DO |
|--|--|
| Battery exhausted. | Replace the battery. |
| Test strip inserted upside down or incompletely. | Insert the test strip with contact bars end first and facing up. |
| Defective meter or test strips. | Please contact customer services. |

2. If the test does not start after applying the sample:

| POSSIBLE CAUSE | WHAT TO DO |
|---|--|
| Insufficient blood sample. | Repeat the test using a new test strip with larger volume of blood sample. |
| Defective test strip. | Repeat the test with a new test strip. |
| Sample applied after automatic switch-off (2 minutes after last user action). | Repeat the test with a new test strip. Apply sample only when flashing "●" appears on the display. |
| Defective meter. | Please contact customer services. |

3. If the control solution testing result is out of range:

| POSSIBLE CAUSE | WHAT TO DO |
|--|--|
| Error in performing the test. | Read instructions thoroughly and repeat the test again. |
| Control solution vial was poorly shaken. | Shake the control solution vigorously and repeat the test again. |
| Expired or contaminated control solution. | Check the expiry date of the control solution. |
| Control solution that is too warm or too cold. | Control solution, meter, and test strips should be at room temperature (20°C to 25°C / 68°F to 77°F) before testing. |
| Defective test strip. | Repeat the test with a new test strip. |
| Meter malfunction. | Please contact customer services. |

SYMBOL INFORMATION

| SYMBOL | REFERENT |
|---|---|
|  | For in vitro diagnostic use only (for use outside the body) |
|  | Do not reuse |
|  | Read instructions before use |
|  | Keep away from sunlight |
|  | Keep dry |
|  | Temperature limitation |
|  | Use by / Expiry date |
|  | Batch code |
|  | Manufacture date |
|  | Manufacturer |
| SN | Serial number |
|  | Caution, consult accompanying documents |
|  | Do not use if package is damaged |
|  | Use within 3 months after opening |
|  | Sterilized using irradiation |
|  | CE mark |

SPECIFICATIONS

Model No.: Dr. Frei MINI

Dimension & Weight: 86,9(L) x 27(W) x 11,3(H) mm, 19g

Power Source: one lithium battery CR2032 3V

Display: LCD

Memory: 20 measurement results with respective date and time

Auto electrode insertion detection

Auto sample loading detection

Auto reaction time count-down

Auto switch-off after 3 minutes without action

Temperature Warning

Operating Conditions:

10°C to 40°C, below 85% R.H. (non-condensing)

Storage / Transportation Condition:

4°C to 40°C, below 95% R.H.

Measurement Units: mg/dL

Measurement Range: 20 to 600mg/dL (1.1 to 33.3mmol/L)

This device has been tested to meet the electrical and safety requirements of: IEC/EN 61010-1, IEC/EN 61010-2-101, IEC/EN 61326-1, IEC/EN 61326-2-6.

*The specifications and the outer design are subject to change for improvement reasons without previous notice.

2 years of guarantee. 10 years of free service.

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■ FORA CARE SUISSE AG, Neugasse 55, 9000 St. Gallen, Switzerland, tel.: 41 71 220 1001,
Manufactured by the order of: Medhouse Swiss GmbH, Tramstrasse 16, 9442 Berneck, Switzerland,